

BESDATA®

WE HAVE BETTER SOLUTIONS

VIDEO LARYNGOSCOPE

VIDEO STYLET

LIGHTED STYLET

REUSABLE VIDEO LARYNGOSCOPE

Stainless Steel Blades



Blade parameter size

Item Number: BD-DF

Blade	Mil0	MAC1	MAC2	MAC3	MAC4	D blade
Application	Neonate	Pediatric	Small Adult	Adult	Large Adult	Difficult airway
Intubation Inner Diameter	2.5-4.0mm	4.0-6.0mm	5.0-6.0mm	6.0-9.0mm	7.5-10.0mm	7.5-10.0mm
Length (cm)	8.2	10.45	12.11	14.25	16.08	13.27
Sterilize Method	2.0%Glutaral fully soak for 10 minutes Or 0.5%-0.6% Cidex OPA for more than 5 minutes					

DISPOSABLE VIDEO LARYNGOSCOPE

Anti-Fog medical-grade polycarbonate



Blade parameter size

Item Number: BD-DP

Blade	Mil0	MAC1	MAC3	MAC4
Application	Neonate	Pediatric	Adult	Large Adult
Intubation Inner Diameter	2.5-4.0mm	4.0-6.0mm	6.0-9.0mm	7.5-10.0mm
Length (cm)	8.31	10.15	13.69	16.08
Sterilize Method	2.0%Glutaral fully soak for 10 minutes Or 0.5%-0.6% Cidex OPA for more than 5 minutes.			

Technical Specification

Parts	Description	Technical Index
Machine	Display	3inch LCD
	Resolution Ratio	640*480 (RGB)
	Viewing angle	≥60°
	Power	<2W
	Display front and back rotation Angle	0°~130 °
	Weight	225g
Camera	Illuminance	≥150 LUX
	Resolution Ratio	1280*720 px
Image/Video function	Image/ Video function	Yes
	Output	USB output, Easy to storage and establish image/video
	Memory type	Micro SD card 8GB
Battery	Battery Type	Rechargeable lithium battery
	Capacity	1350mAh
	Cycle life of battery	>500 times
	Battery working time	>240min
	Charging time	<2(hrs)
Power Adapter	Charging Port	Micro USB
	Input	100~250V, 50 Hz.
	Output	5V1A
Working Condition	Temperature	-5 C ~+50 C
	Humidity	10%~90%
	Air Pressure	860~1060hpa
Transport / Storage Condition	Temperature	-10 C ~50 C
	Humidity	≤93%
	Air Pressure	500~1060hpa


BESDATA® VIDEO Laryngoscope make it better prepared for unexpected tracheal intubations.

- Confidence — MAC-style video laryngoscope blades for a wide range of clinical preferences
- Simplicity — One-button operation, Automatic power off
- Convenience — Anti-Fog Camera, 3.0inch screen reduces blind spot ,250 minutes battery with screen low battery alarming notice

Compared to traditional MAC-style laryngoscopy Blades, Besdata@Video Laryngoscope articulated tip have a bigger angle that allows to raise the epiglottis equal to or greater than 60 °.

LIGHTED INTUBATION STYLET

- Malleable Stylet
- Red-Light Directive
- Medical-Grade Alloy



MALLEABILITY

Success Rate

▼

99%

Time to Intubate

▼

Less than 15 seconds



Item number	Description	O.D. Size	Tube Sizes I.D.
BD-SL-N	Neonate	2mm	2.4-4.0mm
BD-SL-P	Pediatric	3mm	4.0-6.0mm
BD-SL-A	Adult	4mm	4.5-10mm
Sterilize Method	2.0%Glutaral fully soak for 10 minutes Or 0.5%-0.6% Cidex OPA for more than 5 minutes.		


VIDEO LIGHTED INTUBATION STYLET



- Malleable Stylet
- Red-Light Directive
- Medical-Grade Alloy
- Difficult Visualization



Easy & Affordable



MALLEABILITY

The future choice for difficult and routine airways

Item number	Description	O.D. Size	Tube Sizes I.D.
BD-VLS-P	Pediatric	3.9mm	4.5-6mm
BD-VLS-A	Adult	5.5mm	6.0-10mm
Sterilize Method	2.0%Glutaral fully soak for 10 minutes Or 0.5%-0.6% Cidex OPA for more than 5 minutes.		



BESDATA TECHNOLOGY CO.,LTD

Add:2 / F, Building C, Hongyi Industrial Park,
No.4010 Banxuegang Road, Shenzhen, 518129,China
Tel: +86-755-84416070 Fax: +86-755-23039200
support@besdatatech.com www.besdatatech.com