



Your View, We Build.

ViSU LARYNGOSCOPE



Semi-channel


Single-use

One size


Universal system



Single-use Video laryngoscope



SEMI-CHANNEL TECHNOLOGY



High Success rate at 1st-attempt

No need stylet or bougi.

No need to put tube before intubation

Classical SOP


Patented form

SEMI-CHANNEL SOP

	Classical SOP	Channel SOP	Semi-Channel SOP
Need stylet	Yes	No	No
Tube preload	No	Yes	No

APPLICABILITY: ONE SIZE

One device for all adults, which covers **90%** cases. Any ETT 5.5 - 8.5mm.



FULLY DISPOSABLE SCOPE

NO risk of cross-infection
NO cleaning
NO maintenance
NO Batteries

Total single use

Not reusable



Without cower



UNIVERSAL SYSTEM

Fully disposable Video laryngoscope uses USB to connect with Android, Windows or Mac devices.

You can take a picture, record video or share your screen by Wi-Fi to any other phone, pad, laptop or android TV.

PLAG AND PLAY





Photo

Video

Gallery

Settings



Your View, We Build.

Zhuhai Mindhao Medical Technology Co., Ltd.

Address: 4/F, Phase 2, No.2 Chuangxin Si Road, Tangjiawan Town, Hi-tech Zone, Zhuhai 519085, P.R. China

TEL : 0756-3635312

Web: www.mindhao.com